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Corrosion Resistant Intelligent Graphite Digester

Corrosion resistant graphite digester is a kind of pretreatment equipment suitable for sample digestion in high concentration acid mist environment. It can support the simultaneous processing of up to 60 samples, fully meeting the needs of batch sample processing. Step temperature control and automatic sample lifting and lowering function can automatically complete digestion, acid driving, cooling and other work steps, which can minimize the opportunities for laboratory personnel to contact acid gas.

Main features of the product:

○Multiple protection -- corrosion protection and durability

1. Glass ceramic is used as the upper panel of the instrument, which is easy to clean while antiseptic;
2. Advanced high temperature resistant anti-corrosion coating is adopted on the surface of graphite body to effectively prevent strong acid corrosion;
3. The instrument box is made of stainless steel, and the surface is sprayed with anti-corrosion coating, which has excellent corrosion resistance;
4. The electrical system adopts independent sealed space design to better protect electrical components from acid gas corrosion;
5. The instrument adopts wireless communication control mode, which avoids the disadvantages of traditional wired control mode port easy to corrode.

○Wireless control -- safe and convenient

The intelligent graphite digestion instrument supports single machine operation mode and wireless control mode.

The app control software based on Android system can remotely control the device through the mobile terminal. It can easily complete the digestion process without opening the window of the fume hood, reduce the contact between the laboratory personnel and high temperature acid gas, and better protect the health of the experimental personnel.

1. The mobile phone displays the real-time temperature curve to grasp the digestion process more intuitively;
2. The set temperature and holding time of digestion program can be edited and set by mobile phone;
3. It can control the instrument to start and stop digestion;
4. The built-in expert program library can call expert program at any time, create new program, edit and save each program in the library, and call all programs in the library at any time;
5. Beep alarm at the end of program running.

○Automatic lifting -- Sample protection

1. The sample rises automatically to dissipate heat

After the instrument completes the heating program, the sample holder will automatically rise, and the sample will enter into the heat dissipation state, so as to avoid the sample being damaged by the residual temperature of graphite body if the tester fails to take out the sample in time.

2. Sample hovered to near dry

When the digestion work is nearly completed and the sample quantity is small, the bottom of the sample tube can be lifted to the position where it is easy to observe through the lifting function until the sample is almost dry to avoid sample damage.

3. Observe the sample through the window

In the process of digestion, it is often necessary to check the digestion status of the sample. The sample rack can be lifted by wire control switch or app software. Users can observe the sample through the fume hood glass.

○ Uniform temperature heating synchronous digestion

The intelligent graphite digestion instrument uses a unique carbon fiber homogenization heater as the heat source. The sample is heated by isostatic pressing graphite body. The heating speed is faster, the temperature between digestion holes is more uniform, and the service life is longer.

Series	Appearance	Model	Digestion hole specification	Number of digestion holes	temperature range	Temperature difference between holes	Heating level	Power (W)	Dimension (width × depth × height) mm
mini Heater		miniHeater 1530	DiaΦ32mm Hole depth 65mm	15	RT ~ 240°C	±1.5°C (@200°C)	20	1500	290×300×230
		miniHeater 1242	DiaΦ44mm Hole depth 65mm	12	RT ~ 240°C	±1.5°C (@200°C)	20	1500	290×300×230
EH 系列		DXEH-2430	DiaΦ32mm Hole depth 65mm	24	RT ~ 240°C	±1.5°C (@200°C)	20	1800	350×320×230
		DXEH-1542	DiaΦ44mm Hole depth 65mm	15	RT ~ 240°C	±1.5°C (@200°C)	20	1800	350×320×230
		DXEH-3630	DiaΦ32mm Hole depth	36	RT ~ 240°C	±1.5°C (@200°C)	20	2500	450×420×300

			65mm						
		DXEHS-2542	DiaΦ44mm Hole depth 65mm	25	RT ~ 240°C	±1.5°C (@200°C)	20	2500	450×420× 300
		DXEHS-6030	DiaΦ32mm Hole depth 65mm	60	RT ~ 240°C	±1.5°C (@200°C)	20	3500	640×420× 300
		DXEHS-4042	DiaΦ44mm Hole depth 65mm	40	RT ~ 240°C	±1.5°C (@200°C)	20	3500	640×420× 300
		DXEHS-3630	DiaΦ32mm Hole depth 65mm	36	RT ~ 240°C	±1.5°C (@200°C)	20	2500	530×420× 310
		DXEHS-2542	DiaΦ44mm Hole depth 65mm	25	RT ~ 240°C	±1.5°C (@200°C)	20	2500	530×420× 310
		DXEHS-6030	DiaΦ32mm Hole depth 80mm	60	RT ~ 240°C	±1.5°C (@200°C)	20	3500	720×420× 310
		DXEHS-4042	DiaΦ44mm Hole depth 80mm	40	RT ~ 240°C	±1.5°C (@200°C)	20	3500	720×420× 310